


## Bite problems in the morning?

Your muscles have taken another orientation. This might last for a while, but will have no consequences.

Bite for a couple of minutes as shown onto this bite-help strip and the reorientation will be considerably accelerated.



Before and after use, clean bite-help strip with soap, rinse well and store in a dry place.



bite help strip: order number:, 59 60 20, 20 pieces

 **ERKODENT**<sup>®</sup>

Siemensstraße 3 • 72285 Pfalzgrafenweiler  
Tel.: + 49 (0) 74 45/85 01-0 • [info@erkodent.com](mailto:info@erkodent.com)